



**Your Dependable  
Liquid Handling Partner**  
EzMate Automated Pipetting System

# Automated Pipetting System

## EzMate™ 401 / 601

### Your Dependable PCR / qPCR Setup Partner

EzMate™ 401 / 601 is an automated, high-precision pipetting system specifically designed for low-volume PCR / qPCR sample preparation. Lab technicians understand the many challenges associated with setting up PCR and qPCR reactions. Small volumes of reagents need to be dispensed accurately and consistently for acceptable PCR results, and the manual liquid handling process is tedious and time consuming, with many chances for mistakes.

The EzMate™ 401 / 601 was designed as a simple instrument to replace the manual process of PCR / qPCR sample prep. Accuracy, precision and consistency are guaranteed, and reagents will not be wasted on mistakes. Unlike complex, multi-purpose robotic systems, our instrument was designed for the researcher without prior robotics experience. Intuitive set up and programming will have your lab saving time and money almost immediately. Choosing EzMate™ 401 / 601 is a choice of "Working Smart".



#### ► Easy and Convenient to Use

- PC Software : EzStarter™ can be mastered in one hour. No trained technician is required.
- Built-in PCR / qPCR setup protocols can be quickly modified and transferred via USB memory stick.
- 1 and 8-channel, 50µl or 200µl, Automated Pipetting Modules (APMs) can be exchanged without tools.
- 4 or 6 standard microtiter plate / tip rack areas and 2 reagent areas.
- Wide variety of adapters for different types of labware.
- Import pipetting commands via .csv / .txt format files.
- Other Applications : Sequencing setup, HLA typing, SNP detection, DNA normalization

#### ► Easy to Service

- The Automated Pipetting Module (APM) is simple to remove and can be mailed-in for service and calibration.
- PC software is upgradable via the internet.
- The entire system is compact and light-weighted.


#### ► Easy to Afford

- The most affordable Automated Pipetting System available.
- EzTip™ robotic tips are compatible with Beckman's® Biomek® 3000 system.
- The passive cooling CoolBlock™ adapters keep sensitive reagents / samples at 7°C for more than 60 minutes.
- Save reagent costs by reducing human errors and using more dense plates.

#### ► Accurate and Precise


- Each Automated Pipetting Module (APM) is calibrated following ISO-8655 standards.
- Excellent results for qPCR standard curve and replicates (Figure 1).
- Better Precision of pipetting results than manual pipetting (Figure 2).

## Automated Pipetting System



**Figure 1. Excellent Serial Dilution Results**

7ml NIH3T3 Cell cDNA sample is diluted in 21ml water at 1:4 ratio for 4 times. Roche® LightCycler® 480 real-time PCR thermal cycler and Finnzyme DyNAmo® Flash SYBR® Green qPCR kit (F-415L) were used.



**Figure 2. Better Precision Results than Manual Pipetting**

Target: Human GAPDH (top curve) was amplified in replicates of 18µl (master-mix) and 2µl (cDNA) were pipetting into 20µl reaction volume. Roche® LightCycler® 480 real-time PCR thermal cycler and Invitrogen™ Platinum® Taq DNA Polymerase were used.

### 1. Functional Safety Housing

- Safety door detection for emergency stop
- Light-weighted Acrylic material

### 2. Automated Pipetting Module (APM)

- Easy to exchange
- 1/8 channel, 50ml / 200ml models
- Automatic channel and volume identification

### 3. Disposable Used Tip Tray

- Capacity > 400 tips

### 4. 2x Reagent Areas

- For 1.5ml or 2ml micro-centrifuge tubes, 8ml reagent bottle, 15ml or 80ml reagent reservoirs and 8-strip PCR tubes
- CoolBlock™ available for 1.5ml, 2ml micro-centrifuge tubes, 8ml reagent bottle

### 5. 4 / 6-position SBS Worktable

- 1 designated tip rack area
- 2 designated SBS format labware areas
- 1 (3 for 601 systems) Common areas for both labwares and tip racks.
- CoolBlock™ available for 96 / 384-well PCR plates and 20 well adapter for 1.5ml or 2ml micro-centrifuge tubes

### 6. Compact Design

- Small footprint
- Extremely light-weighted

### 7. 15.6" Notebook Computer with EzStarter™ Software

- PCR / qPCR-oriented software
- Pre-run simulation
- Operating System: Windows® XP / 7 / 8 (32 / 64 bit)



### Performance of Automated Pipetting Module (APM)

1 / 8 Channel-Volume: 50µl		
	1µl	50µl
Accuracy (Rel.)	± 7%	± 1%
Precision (Rel. CV)	≤ 7.5%	≤ 0.4%

1 / 8 Channel-Volume: 200µl		
	10µl	200µl
Accuracy (Rel.)	± 3%	± 0.8%
Precision (Rel. CV)	≤ 1%	≤ 0.15%

Note: According to ISO-8655 standards (Gravimetric method), APM is calibrated in temperature (21~25°C, ± 0.5°C) and humidity (60~90%) controlled environment. Twice-distilled water, robotic tips and microbalance were used.

## EzMate™ 401s / 601s

### For High-contamination and Temperature-sensitive Applications


The all new EzMate™ 401s and 601s use the same platform and software of EzMate™ 401 / 601, but have more advanced features added, such as UV lamp, HEPA filters and Active Cooling / Heating modules (ACHM). For high-contamination clinical labs and RNA experiments, the UV lamp and HEPA filter provide a clean working environment. For temperature-sensitive samples and reagents, the Active Cooling / Heating module (ACHM) prevent them from degrading.



### Exceptional Safety

The advanced design of the EzMate™ 401s / 601s delivers improved process safety, giving you peace of mind.

The standard UV and HEPA modules ensure effective worktable decontamination and provide positive clean air pressure throughout PCR setup. The UV safety shield provides protection to human eyes from UV light.



## Automated Pipetting System

To heat or to cool, it's up to you!

With optional Active Cooling / Heating Module (ACHM), you can keep your precious samples and reagents cool and fresh during the experiment setup. You can also heat up your experiment mixture to speed up the reaction. Interchangeable cooling / heating adapter allows you to adopt different types of labware.

The Active Cooling / Heating Module (ACHM) includes the Active Cooling / Heating Control Box, which provides the user an interface, a power supply, and the Active Cooling / Heating Adapters.

### Active Cooling / Heating Module (ACHM)



Active Cooling / Heating Positions: R1 / A / C Areas



384-well PCR plate Active Cooling / Heating Adapter



96-well PCR plate Active Cooling / Heating Adapter



R1 Active Cooling / Heating Adapter



R2 Active Cooling / Heating Adapter



ACHM Control Box

### Technical Specification


UV Lamp	
Irradiation period per run	15 min
Irradiation capacity	> 40 mW/cm <sup>2</sup>
Wavelength	254 nm (UV-C)
HEPA System	
Filter	3M® Air-Mate™ HEPA Filter x 2
Volume flow	> 200 L/min
Active Cooling / Heating Module (ACHM)	
Setting Range	4 ~ 70°C
Max. Heating Rate	5°C/min at RT=22°C
Max. Cooling Rate	2°C/min at RT=22°C
Available Position	R1 / A / C Areas

## Automated Pipetting System

## EzApp™ DNA/RNA Normalization


Having trouble with sample concentration normalization prior to reverse transcription? Want to make sure that all the samples are loaded equally to reduce internal error? EzApp™ DNA/RNA Normalization is here for you. No more tedious calculation. No more pipetting error. All you need to do is providing the original concentration of your DNA/RNA samples, the final concentration, and the volume you need, and EzApp™ will do the rest for you.

## Sample Normalization Procedure




## DNA Normalization Data

## 20 Different Arabidopsis DNA Samples Normalization Data



DNA Normalization  
by EzApp™



## More EzApps are Coming Soon!

We can provide custom-made EzApp™ upon request.



## Product Specification



	401	601	401s	601s
Worktable Capacity				
SBS Format Area	2 ~ 3	2 ~ 5	2 ~ 3	2 ~ 5
Tip Rack Area	1 ~ 2	1 ~ 4	1 ~ 2	1 ~ 4
Reagent Area	2	2	2	2
Function				
Liquid (Sample/Reagent) Transfer (LT)	v	v	v	v
Multiple Dispense (MD)	v	v	v	v
Serial Dilution (SD)	v	v	v	v
Hold (Pause)	v	v	v	v
Mixing (MIX)	v	v	v	v
Loop	v	v	v	v
UV / HEPA	N/A	N/A	v	v
Active Cooling / Heating Module (ACHM)	N/A	N/A	Optional, Up to 3 a reas	Optional, Up to 3 a reas
Automated Pipetting Module	Interchangeable 1 / 8-channel, Maximum volume 50µl / 200µl			
Connection	RS-232, USB 3.0 / 2.0			
Power Supply	100 ~ 240V, 50/60 Hz, 100W			
Size (W x D x H)	590 x 450 x 470 mm	740 x 450 x 470 mm	660 x 450 x 490 mm	810 x 450 x 490 mm
Weight (N.W.)	25 Kg	29 Kg	31 Kg	35 Kg
ACHM Control Box Size (W x D x H)	N/A	N/A	240 x 300 x 230 mm	240 x 300 x 230 mm
ACHM Control Box Weight (N.W.)	N/A	N/A	2 Kg	2 Kg

## Application Note Highlights

- Avoid Tedious Pipetting and Offer Reliability during Antigen Transfer by using EzMate™ 401
- Analyze Copy Number of Mitochondrial DNA by qPCR with the EzMate™ 401
- Clinical Laboratory adopting EzMate™ 401 in Hepatitis B Virus (HBV) Real-time PCR Test
- Tremendous Reagent and Labor Saving  
Saves 40% cost for Real-time PCR setup in high throughput SNP analysis using EzMate™ 401 system and 384-well PCR microplate
- Bone Marrow Registry SBT HLA Typing with EzMate™ 401 Automated Pipetting System
- A Timer Saver in Bone Marrow Transplantation
- Detection of HPV Genotyping by TS-PCR and Biochip with EzMate™ 401
- Improved PCR Uniformity and Sensitivity by Using EzMate™ 401
- Evaluate the Performance of Real-time PCR Instruments by Using EzMate™ 401 Automated Pipetting System
- Plant Leafs Genomic DNA Extraction for Genomic Mapping using EzMate™ 401



## Ordering Information

Catalog No.	Description
<b>EzMate™ 401</b>	
AREZ401-10000	4-position Automated Pipetting System, 15.6" Notebook Computer and EzStarter™ control software
<b>EzMate™ 601</b>	
AREZ601-10000	6-position Automated Pipetting System, 15.6" Notebook Computer and EzStarter™ control software
<b>EzMate™ 401s</b>	
EZ4SU00-100000	4-position Automated Pipetting System with UV / HEPA, 15.6" Notebook Computer and EzStarter™ control software
EZ4SUA0-100000	4-position Automated Pipetting System with UV / HEPA and 2-position Active Cooling / Heating module, 15.6" Notebook Computer and EzStarter™ control software
EZ4SUAC-100000	4-position Automated Pipetting System with UV / HEPA and 3-position Active Cooling / Heating module, 15.6" Notebook Computer and EzStarter™ control software
<b>EzMate™ 601s</b>	
EZ6SU00-100000	6-position Automated Pipetting System with UV / HEPA, 15.6" Notebook Computer and EzStarter™ control software
EZ6SUA0-100000	6-position Automated Pipetting System with UV / HEPA and 2-position Active Cooling / Heating module, 15.6" Notebook Computer and EzStarter™ control software
EZ6SUAC-100000	6-position Automated Pipetting System with UV / HEPA and 3-position Active Cooling / Heating module, 15.6" Notebook Computer and EzStarter™ control software
<b>Interchangeable Automated Pipetting Module</b>	
275-ezar01-01	1-channel, 50ml Pipetting Module (New Bracket)
275-ezar02-01	8-channel, 50ml Pipetting Module (New Bracket)
275-ezar03-01	1-channel, 200ml Pipetting Module (New Bracket)
275-ezar04-01	8-channel, 200ml Pipetting Module (New Bracket)
275-ezar04-11	8-channel, 200ml Pipetting Module for Axygen FX-250-R tips (New Bracket)
<b>Robotic Tips</b>	
T00-ezar00-00	EzTip™ 50µl, non-sterile, 96 tips/rack; 10 racks/pack
T00-ezar01-00	EzTip™ 200µl, non-sterile, 96 tips/rack; 10 racks/pack
T00-ezar02-00	EzTip™ 50µl, sterile, 96 tips/rack; 10 racks/pack
T00-ezar03-00	EzTip™ 200µl, sterile, 96 tips/rack; 10 racks/pack
<b>Accessories</b>	
275-ezar21-00	96 tips rack adapter
275-ezar22-00	Elevated 96-well PCR plate adapter
275-ezar23-00	Elevated 384-well PCR plate adapter
275-ezar24-00	Deep well plate adapter
275-ezar26-00	4 x 2 1.5ml tubes adapter
275-ezar27-00	3 x 2 2ml storage tubes & 1 x 5ml bottle adapter
275-ezar30-00	CoolBlock™ 96 adapter for 96-well PCR plates
275-ezar31-00	CoolBlock™ 384 adapter for 384-well PCR plates
275-ezar34-00	CoolBlock™ R10 adapter for 4 x 2 1.5ml tubes
275-ezar35-00	CoolBlock™ R20 adapter for 3 x 2 2ml storage tubes and 1 x 5ml bottle
275-ezar37-00	CoolBlock™ 20 adapter for 1.5 / 2ml tubes
275-ezar40-00	20-well 1.5ml tubes adapter
275-ezar47-00	80ml reservoir adapter
275-ezar49-00	Disposable 80ml reservoir x 20 pcs
275-ezar51-00	Disposable used tip tray x 10 sets
275-ezar66-00	HLA typing adapter for 60/72-well
275-ezar67-00	HLA typing adapter for 96-well
275-ezar68-00	0.5ml tube adapter for 20-well adapter
275-ezar72-00	3 x 15ml reservoirs adapter
275-ezar73-00	Disposable 15ml reservoir x 9 pcs
275-ezar74-00	3 x 8-strip tubes adapter
275-ezar75-00	96-well conical bottom adapter
275-ezar76-00	Elevated ELISA / cell culture adapter
275-ezar78-00	Illumina 48-well plate adapter
275-ezar79-00	R1 Active Cooling / Heating Adapter
275-ezar80-00	R2 Active Cooling / Heating Adapter
275-ezar81-00	96-well PCR plate Active Cooling / Heating Adapter
275-ezar82-00	384-well PCR plate Active Cooling / Heating Adapter

Trademarks: Beckman®, Biomek® 3000 (Beckman Coulter Inc.); Roche®, LightCycler® 480 (Roche Group); DyNAmo™ (Finnzyme Oy); SYBR® (Molecular Probes Inc.); Invitrogen™ Platinum® (Invitrogen Corp.); Axygen® (Corning Inc.); Illumina® (Illumina, Inc.)  
 Note: Specifications are subjected to change without notice.

## Accessories



1-channel APM



HLA typing adapter for 60 / 72-well



CoolBlock™ 96



8-channel APM



HLA typing adapter for 96-well



CoolBlock™ 384



96 tips rack adapter



Illumina 48-well plate adapter



CoolBlock™ 20



Elevated 96-well PCR plate adapter



4 x 2 1.5ml tubes adapter



CoolBlock™ R10



96-well conical bottom adapter



3 x 2 2ml storage tubes and 1 x 5ml bottle adapter



CoolBlock™ R20



Elevated 384-well PCR plate adapter



3 x 8-strip tubes adapter



96-well PCR plate Active Cooling / Heating Adapter



Elevated ELISA / cell culture plate adapter



80ml reservoir adapter



384-well PCR plate Active Cooling / Heating Adapter



Deep well plate adapter



Disposable 80ml reservoir



R1 Active Cooling / Heating Adapter



20-well 1.5ml tubes adapter



3 x 15ml reservoirs adapter



R2 Active Cooling / Heating Adapter



0.5ml tube adapter for 20-well adapter



Disposable 15ml reservoir



Disposable used tip tray set



ISO 13485