

Centrifuges




The compact and robust low-speed model, 726 406 is well suited for routine work with small sample sizes less than 15 ml tubes. It has an open style of 15 ml angle rotor that can spin at max. 4,000 rpm (2,075 x g).

The product is basically supplied together with an angle rotor and six 15 ml aluminum sleeves, which can be spun immediately upon receipt. Optional tube adaptors for 3 ml, 4 ml and 5 ml tubes are available for easy separation of the small tubes including vacutainers.

It features maximized simplicity in operation and highest safety in everyday work.

It is commonly used in bio-banks, clinical and commercial laboratories for the fractionation of blood samples and simple precipitation of clinical and biological samples.

► Control Panel



Centrifuges

► Features :

- Commonly used in the fractionation of blood samples and simple precipitation of biological samples
- Automatic lid release and vivid end alarm for easy retrieval of samples
- Supplied with an angle rotor including 6 of aluminum sleeves for 15 ml tubes
- CE MDD-certified version, MD406, is also available (if required) **CE 0120**

► Technical Specifications :

Model No.	726 406
Max. RPM	4,000
Max. RCF	2,075 xg
Max. Capacity	6 x 15 ml
Time control	Timed < 100 min or continuous
RPM/RCF conversion	Yes
Noise level (dB)	~ 56
ACC/DEC (sec)	≤ 20 / ≤ 20
Program memory	10
Parameters on display window	RPM (RCF), Min (Sec)
Imbalance cut-off	Yes
Safety lid lock	Yes
Lid drop protection	Yes
Automatic lid release at completion	Yes
Chamber material	Aluminum/Anodizing
Power supply	230 V, 50 Hz (AC 220-230 V, 50/60 Hz; 110 V optional)
Power requirement (VA)	140
Dimension (W x D x H)	296 x 412 x 206 mm
Weight without rotor (kg)	17.5
CE mark	Yes



► Rotors & Adaptors



Aluminum 15ml Sleeve GLB-15A-FA
 Sleeve bore (ΦxL) : 18 x 90
 Radius (mm) : 116
 Max.RPM : 4,000
 Max.RCF (g-force) : 2,075

Angle Rotor, GRA-15-6A

6 x 15 ml
 ∠ 45°
 Hole diameter (mm) : 18
 Max. height for tube it (mm) : 125
 (120 for conical tube)

Tube						
Tube capacity (ml)	3	4	5	10	15	15 ml conical
Adaptor				None	None	None
Cat. No.	GAS-3(15)	GAS-4(15)	GAS-5(15)	-	-	-
Adaptor bore (ΦxL, mm)	13.5 x 61	13.5 x 65	13.5 x 85	-	-	-
Radius (mm)	95.5	98.5	111	115	116	112
Max.RPM	4,000					
Max.RCF (g-force)	1,703	1,762	1,986	2,057	2,075	2,003

Centrifuges




726 416 is a compact bench-top centrifuge with a maximum capacity of 400 ml at a time. It has five versatile rotors, so you can choose the desired rotor, buckets and adaptors according to your application.

With GRS-S-100-4 rotor and universal 100 ml buckets and diverse adaptors, you can enjoy everyday laboratory work to its fullest. The immediate separation of 32 vacutainer samples or small tubes is also feasible with unique vacutainer buckets. An additional offering of microplate rotor extends the 726 416's application.

It is furnished with unique air flow structure and anti-vibration structure. Its ergonomic design with intuitive control panel offers easiness in operation. The SOFT START/STOP function of 726 416 enables gentle acceleration and deceleration, which is effective for vibration-sensitive sample separation such as blood fractionation.

Venchal Scientific 726 416 centrifuge guarantees the highest test standards that you need for your daily work with industry-leading safety and quality.

► Control Panel



Centrifuges

► Features

- Commonly used in the fractionation of blood samples and simple precipitation of biological samples
- A separate button for "Soft" start/stop function to achieve very gentle acceleration and deceleration
- Furnished with unique air flow design and anti-vibration structure to minimize heat and noise generation
- Lid drop protection design allows easiness and safety when loading or unloading samples.
- Automatic lid release and vivid end alarm for easy retrieval of sample
- A large assortment of rotors, buckets and adaptors for microtubes and all sizes of tubes up to 100 ml

- CE MDD-certified version, MD416, is also available (if required)



► Technical Specifications :

Model No.	726 416
Max. RPM	Fixed angle 4,000 , Swing-out 4,000
Max. RCF	Fixed angle 2,700 xg, Swing-out 2,826 xg
Max. Capacity	Fixed angle 16 x 15 ml, Swing-out 4 x 100 ml, 32 vacutainer tube, 4 MTPs
Time control	Timed < 100 min or continuous
RPM/RCF conversion	Yes
Noise level (dB)	~ 56
ACC/DEC (sec)	≤ 20 / ≤ 25
SOFT start/stop	Yes
Program memory	10
Parameters on display window	RPM (RCF), Min (Sec)
Imbalance cut-off	Yes
Safety lid lock	Yes
Lid drop protection	Yes
Automatic lid release at completion	Yes
Chamber material	Aluminum/Anodizing
Power supply	230 V, 50 Hz (AC 220-230 V, 50/60 Hz; 110 V optional)
Power requirement (VA)	340
Dimension (W x D x H)	375 x 480 x 260 mm
Weight without rotor (kg)	23
CE mark	Yes

► Rotors & Adaptors



Swing-out Rotor, GRS-S-100-4

4 loadings

≤ 90°

* GLB-b100-100 includes cap. * (b100) adaptors work fine with cap but (100) adaptors not.



Round 100 ml Bucket with Cap GLB-b100-100

Bucket bore (Φ x L, mm) : 47x99
Radius (mm) : 155.1
Max.RPM:4,000
Max.RCF(g-force) : 2,774
Max.height for tube fit(mm)
: 120(w/ cap) / 130(w/o cap)

Tube										
Tube capacity (ml)	1.5~2.0	2.6~7	4~10	9~15	15 ml conical	15 ml conical	50 ml	50 ml conical	85	100
Tube rack										
Cat. No.	GAM-m2.0-6 (b100)	GAM-7-5 (b100)	GAM-10-5 (b100)	GAM-15-3 (100)	GAM-c15-3 (100)	GAS-c15 (b100)	GAS-50 (b100)	GAS-c50 (b100)	GAS-85 (b100)	GAS-100 (b100)
Rack capacity (ea/4)	6/24	5/20	5/20	3/12	3/12	1/4	1/4	1/4	1/4	1/4
Hole dimension (ΦxL, mm)	11 x 39	13.5 x 60	16 x 60	17.5 x 105	17.2 x 106.5	17.2 x 106.5	29.5 x 95.9	30 x 100	38.5 x 96.4	44.2 x 93

Centrifuges



Round 50 ml Bucket
with Cap

GLB-b50-100

Bucket bore(Φ x L, mm)
: 30.5 x 97
Radius(mm) : 150.9
Max.RPM : 4,000
Max.RCF(g-force) : 2,700
Max.height for tube fit(mm)
: 115(w/ cap) / 125(w/o cap)

Tube				
Tube capacity (ml)	15	15 ml conical	50 ml conical	50
Adaptor				None
Cat. No.	GAS-u15(50)	GAS-c50(50)	-	-
Adaptor bore (Φ x L, mm)	17.3 x 87	29.5 x 17.5	-	-



Dual Round 15 ml Bucket
with Cap

GLB-bd15-100

Bucket bore Φ x(L, mm) : 17 x 98
Radius(mm) : 150.9
Max.RPM : 4,000
Max.RCF(g-force) : 2,700
Max.height for tube fit(mm)
: 123(w/ cap) / 127(w/o cap)

Tube	
Tube capacity (ml)	15
Bucket capacity (ea/4)	2/8



Dual Conical 15 ml Bucket
with Cap

GLB-bdc15-100

Bucket bore Φ x(L, mm)
: 17 x 102.5
Radius(mm): 155.1
Max.RPM: 4,000
Max.RCF(g-force) : 2,774
Max.height for tube fit(mm) : 120

Tube	
Tube capacity (ml)	15 ml conical
Bucket capacity (ea/4)	2/8



Vacutainer 10 ml Bucket

GLB-10-8-100

Bucket bore Φ x(L, mm) : 16 x 89
Radius(mm) : 155.1
Max.RPM : 4,000
Max.RCF(g-force) : 2,774
Max.height for tube fit(mm) : 120

Tube			
Tube capacity (ml)	3	5	10
Bucket capacity (ea/4)	8/32	8/32	8/32
Adaptor			None
Cat No.	GAS-3(10)	GAS-5(10)	-
Adaptor bore (Φ x L, mm)	13 x 57	13 x 82	-
Radius (mm)	121	146	115.1
Max. RPM	4,000		
Max. RCF (g-force)	2,164	2,612	2,774



Swing-out Rotor,
GRS-S-50-4

4 loadings

≤ 90°



Round 50 ml Bucket

GLB-50-50

Bucket bore(Φ x L, mm) : 30.5 x 91
Radius(mm) : 150.9
Max.RPM : 4,000
Max.RCF(g-force) : 2,700
Max.height for tube fit(mm) : 125

Tube				
Tube capacity (ml)	15	15 ml conical	50 ml conical	50
Adaptor				None
Cat. No.	GAS-u15(50)	GAS-c50(50)	-	-
Adaptor bore (Φ x L, mm)	17.3 x 87	29.5 x 17.5	-	-



Conical 50 ml Bucket

GLB-c50-50

Bucket bore Φ x(L, mm)
: 29.5 x 91.5
Radius(mm) : 150.9
Max.RPM : 4,000
Max.RCF(g-force) : 2,700
Max.height for tube fit(mm) : 125

Tube	
Tube capacity (ml)	50 ml conical

Centrifuges



**Round 15 ml Bucket
GLB-15-50**
Bucket bore(Φ x L, mm)
: 17 x 93
Radius(mm) : 150.9
Max.RPM : 4,000
Max.RCF(g-force) : 2,700
Max. height for tube fit(mm)
: 125

Tube



Tube capacity (ml)	15
--------------------	----



**Conical 15 ml Bucket
GLB-c15-50**
Bucket bore(Φ x L, mm)
: 17 x 93
Radius(mm) : 150.9
Max.RPM : 4,000
Max.RCF(g-force) : 2,700
Max.height for tube fit(mm) :
120

Tube



Tube capacity(ml)	15 ml conical
-------------------	---------------



**Dual Round 15 ml Bucket
GLB-d15-50**
Bucket bore (Φx L, mm)
: 17 x 91.5
Radius(mm) : 150.9
Max.RPM : 4,000
Max.RCF(g-force) : 2,700
Max.height for tube fit(mm) :
120

Tube



Tube capacity (ml)	15
Bucket capacity (ea/4)	2/8



**Dual Conical 15 ml Bucket
GLB-dc15-50**
Bucket bore (Φx L, mm)
: 17 x 91.5
Radius(mm) : 150.9
Max.RPM : 4,000
Max.RCF(g-force) : 2,700
Max.height for tube fit(mm) : 120

Tube



Tube capacity (ml)	15 ml conical
Bucket capacity (ea/4)	2/8



**Angle Rotor,
GRA-S-15-16A**
16 x 15 ml
∠ 45°
Hole diameter (mm) : 18
Max. height for tube fit (mm) :
120 for conical tube



**Swing-Out Rotor,
GRS-S-15-6**
6 loadings
∠ 90°



**Aluminum 15ml Sleeve
GLB-15A-FA**
Sleeve bore (Φx L, mm) : 18 x 90
Radius (mm) : 151.5
Max. RPM : 4,000
Max.RCF (g-force) : 2,700

Tube						
Tube capacity (ml)	3	4	5	10	15	15ml conical
Adaptor				None	None	None
Cat. No.	GAS-3(15)	GAS-4(15)	GAS-5(15)	-	-	-
Adaptor bore (ΦxL, mm)	13.5 x 61	13.5 x 65	13.5 x 85	-	-	-
Radius (mm)	134	136	146	150.5	150.9	147.5
Max. RPM	4,000					
Max. RCF (g-force)	2,397	2,433	2,612	2,692	2,700	2,638



**Round 15ml Bucket
GLB-15-15**
Bucket bore(ΦxL):17 x 93.8
Radius(mm):158
Max.RPM:4,000
Max.RCF(g-force):2,826
Max.height for tube fit(mm):125

Tube



Tube capacity (ml)	10	15
--------------------	----	----

**Swing-Out Rotor,
GRS-S-mw-2**

2 loadings
∠ 90°
Area dimension (mm) : 88.5 x 130.3
Max. height for microplate fit (mm) : 35



Plate		
Plate capacity (ml)	MTP	DWP
Holder capacity (ea/2)	1/2	
Radius (mm)	111	
Max. RPM	4,000	
Max. RCF (g-force)	1,986	

Centrifuges




726 6244 is a compact bench-top refrigerated centrifuge with versatility, convenience, stability and sophisticated design and leads to reliable experimental results.

It can spin samples up to a capacity of four 100 ml at 4,200 rpm (3,059 x g) and six 50 ml conical at max. 6,000 rpm (3,960 x g). This model provides three kinds of fixed angle rotors for micro tubes, 15 ml and 50 ml tubes as well as four swing-out rotors. The comprehensive range of rotors and accessories makes it possible for efficient and versatile performances.

The sleek, ergonomic design, vivid display and intuitive control will create a new laboratory environment that is never seen anywhere else. You can choose the color of the control panel between onyx black or clean white.

726 6244 is used in general laboratories and clinical applications including sample separations of cellular materials and blood, and ideal for appropriate experiments for temperature-sensitive samples.


► Control Panel



Centrifuges

► Features

- For up to 6,000 rpm (3,960 xg) with 6 x 50 ml conical tubes and 4,800 rpm (3,888 xg) with 4 x 100 ml tubes
- Commonly used in the fractionation of blood samples and simple precipitation of biological samples
- Fast cooling function to 4°C in 5 minutes for fast startup of cooled samples
- The compressor-off function when lid is open minimizes unnecessary cooling and frosting
- Intuitive control interface with numeric input buttons allows easy setting and editing of running conditions.
- The control panel design can be selected according to user's preference. (The black or white control panel is optionally available.)
- Selectable time mode of "at set SPEED" for counting running time only when reached the set speed
- A separate "PULSE" key for quick spin
- Parameters can be changed during operation.
- Key lock function helps the set parameters secured during operation.
- With motorized lid mechanics, only a gentle touch is sufficient to close and open the lid.
- Lid drop protection design allows easiness and safety when loading or unloading samples.
- A large assortment of rotors, buckets and adaptors for microtubes and all sizes of tubes up to 100 ml



► Technical Specifications :

Model No.	726 6244 W	726 6244 B
Max. RPM	Fixed angle 6,000 RPM, Swing-out 4,800 RPM	
Max. RCF	Fixed angle 3,960 xg, Swing-out 3,888 xg	
Max. Capacity	Fixed angle 6 x 50 ml, 24 x 15 ml, Swing-out 4 x 100 ml, 32 vacutainer tube, 4 MTPs	
Temp. range (°C)	-10°C to +40°C	
FAST cool button	Yes	
Time control	Pulse, Timed < 10 hr or continuous	
Time counting modes	Selectable, at set speed or from starting	
RPM/RCF conversion	Yes	
Noise level (dB)	~ 60	
ACC/DEC	9/10 steps	
Program memory	100	
Parameters on display window	RPM (RCF), Oper Status, Lid Open/Close, Hr:Min, Temp, ACC, DEC	
Display	LCD White Display	LCD Black Display
Automatic rotor identification	Yes	
Imbalance cut-off	Yes	
Safety lid lock	Yes	
Lid drop protection	Yes	
Motorized lid open & close	Yes	
Power supply	230 V, 50 Hz (AC 220-230 V, 50/60 Hz; 110 V optional)	
Power requirement (VA)	900	
Dimension (W x D x H)	584 x 535 x 317 mm	
Weight without rotor (kg)	65	
CE mark	Yes	

Centrifuges



Swing-out Rotor, GRS-SR-50-4

4 loadings
≤ 90°



Round 50 ml Bucket GLB-50-50

Bucket bore(Φ x L, mm)
: 30.5 x 91
Radius(mm) : 150.9
Max.RPM : 4,800
Max.RCF(g-force) : 3,888
Max.height for tube it(mm) :125

Tube				
Tube capacity (ml)	15	15 ml conical	50 ml conical	50
Adaptor			None	
Cat. No.	GAS-μ15(50)	GAS-c50(50)	-	
Adaptor bore (ΦxL, mm)	17.3 x 87	29.5 x 17.5	-	



Conical 50 ml Bucket GLB-c50-50

Bucket bore Φ x L, mm)
: 29.5 x 91.5
Radius(mm) : 150.9
Max.RPM : 4,800
Max.RCF(g-force) : 3,888
Max.height for tube it(mm) : 125

Tube	
Tube capacity (ml)	50 ml conical



Round 15 ml Bucket GLB-15-50

Bucket bore(Φ x L, mm)
: 17 x 93
Radius(mm) : 150.9
Max.RPM : 4,800
Max.RCF(g-force) : 3,888
Max.height for tube it(mm) :125

Tube	
Tube capacity (ml)	15



Conical 15 ml Bucket GLB-c15-50

Bucket bore(Φ x L, mm)
: 17 x 93
Radius(mm) : 150.9
Max.RPM : 4,800
Max.RCF(g-force) : 3,888
Max.height for tube it(mm) :120

Tube	
Tube capacity(ml)	15 ml conical



Dual Round 15 ml Bucket GLB-d15-50

Bucket bore(Φ x L, mm)
:17 x 91.5
Radius(mm) : 150.9
Max.RPM : 4,800
Max.RCF(g-force) : 3,888
Max.height for tube it(mm) :120

Tube	
Tube capacity (ml)	15
Bucket capacity (ea/4)	2/8



Dual Conical 15 ml Bucket GLB-dc15-50

Bucket bore Φ x L, mm)
: 17 x 91.5
Radius(mm) : 150.9
Max.RPM : 4,800
Max.RCF(g-force) : 3,888
Max.height for tube it(mm) : 120

Tube	
Tube capacity (ml)	15 ml conical
Bucket capacity (ea/4)	2/8

Centrifuges



Swing-out Rotor, GRS-SR-100-4

4 loadings
 $\angle 90^\circ$



Dual Conical 15 ml Bucket with Cap GLB-bdc15-100

Bucket bore(Φ x L, mm)
: 17 x 102.5
Radius(mm) : 155.1
Max.RPM : 4,200
Max.RCF(g-force) : 3,059
Max.height for tube it(mm) : 120

Tube



Tube capacity
(ml)

15 ml conical

Bucket capacity
(ea/4)

2/8

* GLB-b100-100 includes cap. * (b100) adaptors work in with cap but (100) adaptors not.

Round 100 ml Bucket with Cap GLB-b100-100 Bucket bore (Φ x L, mm) : 47 x 99 Radius (mm) : 155.1 Max.RPM : 4,200 Max.RCF(g-force) : 3,059 Max.height for tube it(mm) : 120(w/ cap) / 130(w/o cap)	Tube										
	Tube capacity (ml)	1.5~2.0	2.6~7	4~10	9~15	15 ml conical	15 ml conical	50 ml	50 ml conical	85	100
	Tube rack										
	Cat. No.	GAM-m2.0-6 (b100)	GAM-7-5 (b100)	GAM-10-5 (b100)	GAM-15-3 (100)	GAM-c15-3 (100)	GAS-c15 (b100)	GAS-50 (b100)	GAS-c50 (b100)	GAS-85 (b100)	GAS-100 (b100)
	Rack capacity (ea/4)	6/24	5/20	5/20	3/12	3/12	1/4	1/4	1/4	1/4	1/4
	Hole dimension (Φ x L, mm)	11 x 39	13.5 x 60	16 x 60	17.5 x 105	17.2 x 106.5	17.2 x 106.5	29.5 x 95.9	30 x 100	38.5 x 106.4	44.2 x 106.4



Round 50 ml Bucket with Cap GLB-b50-100 Bucket bore (Φ x L, mm) : 30.5 x 97 Radius(mm) : 150.9 Max.RPM : 4,200 Max.RCF(g-force) : 2,976 Max.height for tube it(mm) : 115(w/ cap) / 125(w/o cap)	Tube				
	Tube capacity (ml)	15	15 ml conical	50 ml conical	50
	Adaptor			None	
	Cat. No.	GAS- μ 15(50)	GAS-c50(50)	-	
	Adaptor bore (Φ x L, mm)	17.3 x 87	29.5 x 17.5	-	



Dual Round 15 ml Bucket with Cap GLB-bd15-100

Bucket bore Φ x L, mm) : 17 x 98
Radius(mm) : 150.9
Max.RPM : 4,200
Max.RCF(g-force) : 2,976
Max.height for tube it(mm) : 123(w/ cap) / 127(w/o cap)

Tube



Tube capacity
(ml)

15

Bucket capacity
(ea/4)

2/8

Centrifuges



Vacutainer 10 ml Bucket GLB-10-8-100

Bucket bore $\Phi \times L$ (mm) : 16 x 89
Radius(mm) : 155.1
Max.RPM : 4,200
Max.RCF(g-force) : 3,059
Max.height for tube it(mm) : 120

Tube			
Tube capacity (ml)	3	5	10
Bucket capacity (ea/4)	8/32	8/32	8/32
Adaptor			None
Cat No.	GAS-3(10)	GAS-5(10)	-
Adaptor bore ($\Phi \times L$, mm)	13 x 57	13 x 82	-
Radius (mm)	121	146	155.1
Max. RPM	4,200		
Max. RCF (g-force)	2,386	2,879	3,059



Dual Conical 50 ml Bucket GLB-dc50-100

Bucket bore $\Phi \times L$ (mm) : 30 x 89
Radius(mm) : 155.1
Max.RPM : 4,200
Max.RCF(g-force) : 3,059
Max. height for tube it(mm) : 125

Tube	
Tube capacity (ml)	50ml conical
Bucket capacity (ea/4)	2/8

Swing-Out Rotor, GRS-SR-15-6

6 loadings
 $\angle 90^\circ$



Round 15ml Bucket
GLB-15-15
Bucket bore($\Phi \times L$, mm): 17 x 93.8
Radius(mm):158
Max.RPM:4,400
Max.RCF(g-force) : 3,420
Max.height for tube it(mm) : 125

Tube		
Tube capacity (ml)	10	15



Swing-Out Rotor, GRS-SR-mw-2

2 loadings
 $\angle 90^\circ$
Area dimension (mm) : 88.5 x 130.3
Max. height for microplate it (mm) : 35

Plate		
Plate capacity (ml)	MTP	DWP
Holder capacity (ea/2)	1/2	
Radius (mm)	111	
Max. RPM	4,000	
Max, RCF (g-force)	1,986	

Centrifuges



Fixed Angle Rotor, GRF-SR-50-6

6 x 50 ml
 $\angle 25^\circ$
 Hole diameter (mm) : 29.5
 Max. height for tube it (mm) : 120

Tube				
Tube capacity (ml)	10	15	15 ml conical	50
Adaptor				None
Cat. No.	GAS-15(50)		GAS-c15(50)	-
Adaptor bore ($\Phi \times L$, mm)	17 x 94		17 x 105	-
Radius (mm)	87.9	88.9	91	95
Max.RPM	6,000			
Max.RCF (g-force)	3,538	3,578	3,663	3,824



Fixed Angle Rotor, GRF-SR-c50-6

6 x 50 ml conical
 $\angle 25^\circ$
 Hole diameter (mm) : 29.5
 Max. height for tube it (mm) : 125

Tube		
Tube capacity (ml)	15 ml conical	50 ml conical
Adaptor		None
Cat. No.	GAS-c15(c50)	-
Adaptor bore ($\Phi \times L$, mm)	17 x 105	-
Radius (mm)	93.6	98.4
Max.RPM	6,000	
Max.RCF (g-force)	3,767	3,960



Angle Rotor, GRA-SR-15-24

24 x 15 ml
 $\angle 45^\circ$
 Hole diameter (mm) : 18
 Max. height for tube it (mm) : 125
 (120 for conical tube)

Tube						
Tube capacity (ml)	3	4	5	10	15	15 ml conical
Adaptor				None	None	None
Cat. No.	GAS-3(15)	GAS-4(15)	GAS-5(15)	-	-	-
Adaptor bore ($\Phi \times L$, mm)	13.5 x 61	13.5 x 65	13.5 x 85	-	-	-
Radius (mm)	116.3/132.3	117.7/133.7	136.8/152.8	137.5/154.5	138/155	133/149
Max.RPM	4,000					
Max.RCF (g-force)	2,080/2,367	2,105/2,463	2,447/2,733	2,460/2,764	2,469/2,773	2,379/2,665



Aluminum 15 ml Sleeve **GLB-15A-FA**
 Sleeve bore ($\Phi \times L$, mm) : 18 x 90
 Radius(mm) : 138/155
 Max. RPM : 4,000
 Max.RCF(g-force) : 2,469/2,773



Fixed Angle Rotor, GRF-SR-m2.0-24

24 x 1.5/2.0 ml
 $\angle 45^\circ$
 Hole diameter (mm) : 11.1
 Max. height for tube it (mm) : 52

Tube			
Tube capacity (ml)	0.2	0.5	1.5/2.0
Adaptor			None
Cat. No.	GAS-m0.2(2)	GAS-m0.5(2)	-
Adaptor bore ($\Phi \times L$, mm)	6.5 x 23	8 x 31	-
Radius (mm)	67	73	84.3
Max.RPM	6,000		
Max.RCF (g-force)	2,697	2,938	3,393