

## Venpet

### Single and Multichannel Electronic Pipettes

- Precise Pipetting
- Ergonomic Force-saving Design
- 360° revolving collar
- Reduce Stress
- Fatigue From Repetitive Pipetting



- **Automatic Switch-On Calibration**  
for high accuracy and precision
- **Unique Force-Saving Design**  
to reduce stress and fatigue from repetitive pipetting
- **360° Revolving Collar**  
for greatest positioning comfort
- **9 Memory Settings**  
to save protocol set-up time
- **User-Friendly Interface**  
and easy-to-operate buttons
- **High Capacity Lithium Battery**  
for high stamina continuous use
- **5 Speeds for Aspiration and Dispensing**  
depending on liquid viscosity



6300020



6300200



6301000



6300208



6302008

## ► Specifications & Ordering Information :

Cat. No.	Description	Volume Range	Increment	Accuracy Rel. % ( $\pm$ )	Precision Range Rel. CV % ( $\leq$ )
6300020	1-Channel	2-20 $\mu$ l	0.1 $\mu$ l	1-4	0.4 - 2.5
6300200	1-Channel	10-200 $\mu$ l	1 $\mu$ l	0.6 -2.5	0.15 - 1
6301000	1-Channel	100-1000 $\mu$ l	1 $\mu$ l	0.5-2.5	0.15 - 0.6
6300208	8-Channel	2-20 $\mu$ l	0.1 $\mu$ l	1-4	0.5 - 2.5
6302008	8-Channel	10-200 $\mu$ l	1 $\mu$ l	0.6-3	0.15 - 1

## Accessories

Cat. No.	Description
6300000-B	High-Capacity Lithium Battery
6300000-C	Standard AC-DC Adaptor and Optional Charger
6300000-S	Stand for 3 Units

Specifications are subject to change without prior notice