


Elisa Microplate Washer



Elisa microplate washer, with 8 or 12 washing probes suitable with 96-well and 48-well microplates.

► Features :

- Selective washing parameters for different effects.
- External programming/downloading available through RS-232 ports.
- Regulates rinse position of cleaning head according to user requirements.
- Controlled by microprocessor, precise and easy to use, washer can automatically accomplish the single and multi-row microplate washing.

► Specifications :

Model No.	512-201
Microplate Types	48-well, 96-well
Washing	Automatic washing, 5 ways
Washing the row	1 row per time
No. of Microplates	1- 12 rows
Average Residue	≤ 2ul per well
Dispense Precision	CV≤3%
Shaking	5 ways
Memory	100 programs
Wash Cycles	1-20
Soak Time	Programmable in seconds, up to 20 minutes
Anti-overflow	Excess fluid is automatically discharged
Power Supply	AC110V/220V ±10%, 50/60Hz
Power	≤80W
Packing Dimension (WxDxH)	716 x 696 x300mm
Gross Weight	18 Kg
Accessories	Power Cord, Fuse Tube, Wash Bottle, Waste Bottle, Pass Needle, Silicone Tube, 8 Needles Washing Head, 12 Needles Washing Head, Manual, Warranty Card

► Ordering Information :

Model No.	Description
512-201	Elisa Microplate Washer, Microplates 48-well, 96-well