

## Ice Flaker



### FEATURES :

- High-quality stainless steel is used for their shells, they are anti-corrosive and durable, and the all-in-one independent structure makes them simple and compact and space-saving.
- Compressor/Refrigerant : High Efficient Fluorine-free R134a
- Motor Cooling : Air Cooling
- Kind of Ice : Irregular small granular snowflake crushed ice
- Their boxes are protected with fluorine-free foam thermal cover featuring good insulation effect. Their internal bladders are of fluorine-free bacterial inhibitory-type. They are both energy-saving and environmental friendly.

### SPECIFICATIONS :

Model No.	934-20	934-30	934-50	934-85	934-100	934-130	934-200	934-250	934-300
Ice making capacity(kg/24h)	20	30	50	85	100	130	200	250	300
Ice storage capacity(kg)	10	10	15	25	25	35	55	55	55
Ice shape	Irregular small granular snowflake crushed ice								
Input current(w)	180	180	280	420	420	685	805	1370	1370
Ambient Temp.	10~25°C								
Water Temp.	10°C								
Power Supply	AC220V±10%,50/60HZ; 110V±10%,60HZ								
Refrigeration type	Air Cooling								
Tank shell	Stainless steel								
Compressor / Refrigerant	High efficient and Fluorine free/R134a								
Standard Accessories	Water inlet pipe-1 no., Drainage pipe - 1no, Ice spoon - 1no, Sealing washer - 2 no								
Optional Accessories	Filter (Water pressure needs to be within 1.5~3Kg/cm <sup>2</sup> )								
Net weight(kg)	32	32	40	57	60	62	85	91	100
Gross weight(kg)	36	36	45	62	67	68	91	98	106
External Size(WxDxH)mm	330x503x607		380x543x722	398x543x837	635x611x945	635x611x945	680x641x1102	680x641x1102	680x641x1102
Package Size(WxDxH)mm	400x574x655		457x621x770	475x622x886	712x690x994	712x690x994	757x720x1150	757x720x1150	757x720x1150
Applicable area	Various types of experiments, supermarkets, hotels and other seafood products stores etc.								

### Ordering Information :

Model No.	Description
934-20	Ice Making Capacity 20 kg/24h, Ice Storage Capacity 10 Kg
934-30	Ice Making Capacity 30 kg/24h, Ice Storage Capacity 10 Kg
934-50	Ice Making Capacity 50 kg/24h, Ice Storage Capacity 15 Kg
934-85	Ice Making Capacity 85 kg/24h, Ice Storage Capacity 25 Kg
934-100	Ice Making Capacity 100 kg/24h, Ice Storage Capacity 25 Kg
934-130	Ice Making Capacity 130 kg/24h, Ice Storage Capacity 35 Kg
934-200	Ice Making Capacity 200 kg/24h, Ice Storage Capacity 55 Kg
934-250	Ice Making Capacity 250 kg/24h, Ice Storage Capacity 55 Kg
934-300	Ice Making Capacity 300 kg/24h, Ice Storage Capacity 55 Kg