



Features

- LED 2-screen digital display actual and setting R.P.M.
- Designed for continuous-operation orbital shaking L
- Differential clamps available
- Three-point bearing transmission is adopted suitable for heavy loads and continuous duty
- Rubber non-slip mat to prevent slipping
- Build-in timer is up to 999/h,m,s.and alarm when time out
- Uses brushless AC motor with spring



► Technical Specifications :

Model No.	219500	219505	219501
Panel Type	Digital		
Motion	Orbital		Rocking
Rate of Shaking	20~250 r.p.m	20~150 r.p.m	6~60 r.p.m
Timer	999 h / 999 min / 999 sec		
Shaking Width	25 mm		35 mm
Shaking Plate	330 x 280 mm		
Power Supply	Power Voltage : 110V/220V 50Hz/60Hz		
Outside Dimension (WxDxH)	330 x 280 x 170 mm	330 x 280 x 370 mm	330 x 280 x 250 mm
Weight	11 kg	12 kg	10 kg

At any time, the manufacturer may bring modifications to the model above described, for technical or commercial reasons. The picture does not necessary shows the product in the standard version.