

# Incubator

## Orbital Shaking Incubator - Bench Top

Orbital Shaking Incubator

### Features

- Interior is made of Stainless Steel SUS304 and Exterior is anti corrosion powder coated
- Touch Screen Panel .
- Overheat and over low protector, built-in circulating fan for temperature stability.
- With safety switch to stop shaking when opening door
- Tempered Safety Glass Window for operation view
- Uses brushless AC motor for long durability and stable performance of the shaking system
- Shaking smoothly and quietly

1254504



### Specifications

Model No.	1254004	1254504	1254204
Type	Non- Refrigerated		
Door Type	Flip-up door		Upright - Front
Temp. Range	Ambient +5°C ~70°C		
Shaking Speed	20~300 r.p.m		
Temp. Accuracy	±0.1		
Temp. Control	PID		
Temp. Display	LED		
Shaking Width	25 mm		
Timer	999 h / 999 min / 999 sec		
Heater	700W		
Shaking Plate Size	400 x 300 mm	450 x 450 mm	400 x 300 mm
Volume	66 L	102 L	65 L
Inside Dimension (WxDxH)	455 x 390 x 380 mm	500 x 540 x 380 mm	450 x 390 x 370 mm
Outside Dimension (WxDxH)	460 x 510 x 540 mm	510 x 650 x 540 mm	520 x 550 x 610 mm
Weight	38 kg	50 kg	55 kg
Power Supply	110V/220V, 50Hz/60Hz		

### Capacity for Flasks

Flask Clamp	50ml	125ml	250ml	500ml	1000ml	2000ml
1254004	35	20	12	6	4	2
1254204	35	20	12	6	4	2
1254504	64	36	25	16	9	5

At any time, the manufacturer may bring modifications to the model above described, for technical or commercial reasons. The picture does not necessary shows the product in the standard version.