

Circulating Water Bath



FEATURES :

- Chamber is made of seamless stainless steel and Exterior body is powder coated
- Microprocessor PID Controller
- Overheat Protection
- Water Level Protector
- LED Digital Display
- No Welded Corners
- The water bath with circulator using electric heat tube to unify temperature inside
- Add exterior circulator as require
- The sensor is PT-100 Ω
- Gable lid Included
- Heater of stainless steel
- External Circulating Pump is Available
- Heating circulator sold separately for your option

Specifications :

Model No.	10012	10013
Temp. range	Ambient +5°C~100°C	Ambient +5°C~100°C
Temp. Accuracy	$\pm 0.1^{\circ}\text{C}$	
Temp. Control	Micro processor PID Controller	
Temp. Sensor	PT-100 Ω	
Temp. Display	LED	
Capacity	10 L	20 L
Circulation Pump	YES	YES
Power Supply	220V/50Hz 1Phase, can be changed as 220V/60Hz or 110V 50V/60Hz	
Inside Dimension (WxDxH)	295 x 235 x 150 mm	495 x 295 x 150 mm
Outside Dimension (WxDxH)	340 x 280 x 310 mm	540 x 340 x 315 mm
Weight	7 kg	9 kg

Ordering Information :

Model No.	Description
10012	Circulating Water Bath, Capacity 10 Ltrs
10013	Circulating Water Bath, Capacity 20 Ltrs