

General Purpose Water Bath

Temperature Range Ambient +5°C ~ 99°C



► Features :

- Microprocessor PID Controller ensures the most accurate temperature control
- Temperature Control up to 5°C to 99°C
- Display 4 LED digital display
- Exterior provides durability with high quality stainless steel corrosion resistant bath
- Seamless stainless steel bath
- Over Temperature Protection capability ensures user safety
- Auto-tune enables optimized temperature display
- Timer up to 99 hours 59 minutes
- Error Notification in case of temperature sensor malfunction
- Perforated cover protects heater and sensor from unexpected damage

► Specifications :

Cat. No.	5206	5211	5214	5222
Temp. Range	Ambient 5°C~99°C			
Temp. Accuracy	Set Accuracy $\pm 0.1^{\circ}\text{C}$			
Uniformity	$\pm 1.0^{\circ}\text{C}$			
Capacity	6 Liter	11 Liter	14 Liter	22 Liter
Controller	Microprocessor Digital PID Controller			
Display	LED 4 Digit Display			
Timer	99 Hours 59 Minutes			
Interior	Seamless Stainless Steel (STS 304)			
Exterior	Stainless Steel (STS 430)			
Safety	Over Temperature Cut Off			
Power Supply	110V, 60Hz or 220V, 50/60Hz			
Inner Dimension (WxDxH) mm	300x155x150	300x240x150	325x298x150	500x295x150
Outer Dimension (WxDxH) mm	327x178x280	327x265x280	354x325x280	537x327x280